

SOL-CARE™ Safety Multi-Sample Needle



SOL-CARE™ Safety Multi-Sample Needle with Pre-Attached Holder



The SOL-CARE™ Safety Multi-Sample Blood Collection Needle is simple and intuitive to use, making conversion straightforward

Helps reducing the risk of needlestick injury with a pivoting shielding technology

- 2 ways to activate the cover, so users can adopt their preferred technique
- An audible click confirms that the cover has safely been locked in place
- Thin-wall needle provides effective flow rate to reduce the risk of blood hemolysis
- Low-angle bevel design helps to minimise pain and discomfort for patients
- Compatible with standard holders
- Pre-attached holder version is compatible with commonly used blood evacuation tube systems



SOL-CARE™ Safety Multi-Sample Needle

SOL-CARE™ Safety Multi-Sample Needle with Pre-Attached Holder



SOL-CARE™ Safety Multi-sample Needle

Reference number	Needle outer diameter		Needle length	
	Gauge	Millimeters	Inches	Millimeters
SMSN22125	22G	0.7 mm	1 1/4"	32 mm
SMSN2215	22G	0.7 mm	1 1/2"	38 mm
SMSN21125	21G	0.8 mm	1 1/4"	32 mm
SMSN2115	21G	0.8 mm	1 1/2"	38 mm
SMSN20125	20G	0.9 mm	1 1/4"	32 mm
SMSN2015	20G	0.9 mm	1 1/2"	38 mm

Sales Unit: 48 pcs/Box; 480 pcs/Case

SOL-CARE™ Safety Multi-sample Needle with Pre-attached Holder

Reference number	Needle outer diameter		Needle length	
	Gauge	Millimeters	Inches	Millimeters
SMSNH22125	22G	0.7 mm	1 1/4"	32 mm
SMSNH2215	22G	0.7 mm	1 1/2"	38 mm
SMSNH21125	21G	0.8 mm	1 1/4"	32 mm
SMSNH2115	21G	0.8 mm	1 1/2"	38 mm
SMSNH20125	20G	0.9 mm	1 1/4"	32 mm
SMSNH2015	20G	0.9 mm	1 1/2"	38 mm

Sales Unit: 50 pcs/Box; 400 pcs/Case

Sol-Millennium Europe Sp. z o.o
 Twarda 18, 00-105 Warsaw, Poland
 Tel: +48 22 202 67 23
 info.EU@sol-m.com

sol-m.com

CORP-BC-053-0519

For more information about
**SOL-CARE™ Safety Multi-Sample
 Needle, please contact your local
 SOL-MILLENNIUM® representative.**

